

Diuron plus MCPA

See individual products, Diuron flowable and MCPA amine for further information.

Rates quoted in HerbiGuide are for a 1:1 mixture of Flowable Diuron 500 plus MCPA amine 500. A rate of 500 mL of the mixture is equivalent to 250 mL of Flowable Diuron 500 plus 250 mL of MCPA amine 500. Similarly 500 mL of the mixture is equivalent to approximately 140 grams of Diuron 900 granules plus 500 mL of MCPA amine 500.

This mixture is synergistic so lower rates of the mixture are just as effective as higher rates of the individual components.

Acknowledgments:

Collated by HerbiGuide. Phone 08 98444064 for more information.